
BIARTICULAR

Cabeza Biarticular
Biarticular Head

 **surgival**
GRUPO COSÍAS

Cabeza Biarticular

Biarticular Head

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The **Biarticular Head**, also named **Bipolar or Bicentric**, was conceived fundamentally to reduce the friction that developed, in most cases and in an advanced form, in the cartilage acetabular to the direct contact of the metallic head of an endoprosthesis in a classic hemi arthroplasty.

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The **Biarticular Head** consists of two parts: a metallic **Cup** that is in contact with the acetabulum and one **Insert** in polyethylene that adapts inside the metallic cup whose interior fixes the head of the femoral stem. In this way there are two articulations made, one between acetabulum and cup, and another between the inset an femoral head.

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The **Cup** is made of stainless steel (AISI 316 LVM), and is presented in 15 sizes from 41 to 55 mm, with a cadence of 1 mm, in function of the diameter of the femoral head.

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
The **Insert** is made of polyethylene of a ultra high molecular weight (UHMWPE), that is translated in a minor “debri” by stress. It is presented in four sizes with an internal diameter of 28 mm.

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1 • The way of beginning is to the discretion of the surgeon.

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4 • Opening and exploration of the medular cavity with a narrow and long spoon.

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5 • Introduction of the rasp in order to make the femoral cavity of the femoral prosthesis.

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7 • On the definitive stem we can carry out a trial reduction using the trial cups and trial heads, the instrumentation set includes trial 28 mm diameter heads and taper 12/14.

Note: If the chosen stem is SURGIVAL, we will be able to carry out previously trial reductions on heads and cups even on the rasp.

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8 • Assembling the different components following the next steps:



1st Step: To introduce the head of 28 mm diameter on pressure in the polyethylene insert.

2nd Step: The block head-insert is introduced in the Biarticular cup, check it is correctly fit together.

Note: The assembling of the femoral head, insert and cup is advisable to carry out on the instrumentation table. In the cases that are considered opportune it is also possible to assemble directly on the neck of the femoral stem.

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9 • Impactation of the assembling Biarticular components on the stem and articulation reduction.

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Ref. A1519541	Ø 41 mm
Ref. A1519542	Ø 42 mm
Ref. A1519543	Ø 43 mm
Ref. A1519544	Ø 44 mm
Ref. A1519545	Ø 45 mm
Ref. A1519546	Ø 46 mm
Ref. A1519547	Ø 47 mm
Ref. A1519548	Ø 48 mm
Ref. A1519549	Ø 49 mm
Ref. A1519550	Ø 50 mm
Ref. A1519551	Ø 51 mm
Ref. A1519552	Ø 52 mm
Ref. A1519553	Ø 53 mm
Ref. A1519554	Ø 54 mm
Ref. A1519555	Ø 55 mm

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C/C cono 12/14 Short Neck 12/14 taper	Ref. A1536040
C/M cono 12/14 Medium Neck 12/14 taper	Ref. A1536041
C/L cono 12/14 Long Neck 12/14 taper	Ref. A1536042
C/XL cono 12/14 Extra Long Neck 12/14 taper	Ref. A1536043

3 Discos probadores (par/impar)
Tester disc (even/odd)

Ref. A1519020

4 Posicionador extractor
Positioner-extractor

Ref. A1519021

5 Cincel hueco
Mortising chisel

Ref. A1160020

6 Extractor P/Insertos biarticular
Extractor for biarticular insert

Ref. A1519022

Caja vacía 1 Fil 60 x 30 x 11
Empty container 1 Fil 60 x 30 x 11

Ref. B4500001

Distribuidor oficial de

Tecnología Europea de Vanguardia
Advanced European Technology

Industria del Agave 159-A
Zapopan Industrial Norte
45130, Zapopan, Jalisco
33 2471 0012

spinetech@spinetech.com.mx

www.spinetech.com.mx / www.surgival.com