KAREY®

Vástago KAREY®
KAREY® Stem

surgival
GRUPO COSÍAS

Características técnicas Technical details

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The cementless KAREY- HA primary stem is manufactured in titanium alloy (Ti6Al4V).

Straight stem of quadrangular section divided into two portions:

a) Proximal or metaphyseal, trapezoidal in shape on sagittal and frontal planes, obtaining three-dimensional primary stability and avoiding component subsidence.

b) Distal or diaphyseal, progressively subsized to reduce implant rigidity, improving contact with the cancellous bone, enhancing bone growth and secondary stability of the implant.

The progressive subsizing of the distal portion to avoid diaphyseal filling, reducing the "shielding stress".

The HA coating is applied using the latest technique of Air Plasma Spraying, with a thickness of 130 micron over the whole implant.

Both the HA coated cementless stem, and the KAREY-C cemented stem are implanted using the same instrumentation.

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The metaphyseal semicircular macrostructures increase the contact surface, avoiding component subsidence and having an antivaro effect.

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The longitudinal grooves increase the bone-implant contact surface, avoiding rotation and reducing distal rigidity with respect to the metaphyseal portion.

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Surgical approach • The objective of any surgical hip approach is correct visualization of the acetabulum and proximal femur. The choice of surgical hip approach will be based on the surgeon's preferences

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1 · Femoral neck osteotomy

With the femoral template on the preoperative X-ray, we will determine the proper femoral-neck level of resection. The osteotomy will be performed at a 45° dihedral angle (Fig. 1).

Fig. 2

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Fig. 4

2 · Preparation of the femur

Start opening the femoral canal using a hollow box chisel (Fig. 2), removing enough bone from the medial wall of the greater trochanter parallel to the femoral axis.

Insert the manual conical reamer (Fig. 3), as laterally as possible and parallel to the femoral axis to find the femoral canal.

With the smallest rasp mounted on the rasp-handle (Fig. 4), start the reaming, inserting it as far as the last line of teeth. This operation must be repeated increasing rasp size sequentially until the endomedullar cortical chatter is found. The anteversion is established with these rasping operations.

Fig. 5

Fig. 6

Fig. 7

Fig. 8

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3 · Trial reduction

With the rasp inserted in the femoral canal, carry out the reduction with the trial heads for rasps, in order to select the proper neck length (Fig. 5). Finally, once the neck length has been chosen, stability, tightness and leg length are tested.

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4 · Femoral component insertion

The definitive stem will correspond to the same size as the last rasp used. Mount the stem on the specific insertion handle (Fig. 6), applying impactation until it is inserted to the required depth, at this point the insertion handle is removed. With the final impactor the definitive stem can be impacted in a complementary way, if considered necessary (Fig. 7).

The trial femoral head for stem permit to do a new reduction to assure the length of the neck.

Finally we place the definitive head chosen on the stem and the articulation is reduced with the help of the impactor-reducer (Fig. 8).

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Note: The same instrumentation is used to implant both the KAREY HA stem and the KAREY C stem using the same technic shown. However, take into consideration that for the implantation of the KAREY-C one must consider the steps specific to all cemented stems: cavity cleaning, preparation and insertion of the bone cement, distal obturators, etc.


Instrumental Instruments

Set completo instrumental vástago KAREY
KAREY Stem instrumentation complete set

Ref. F0005005

Bandeja superior
Upper tray

| | |
|---------------|-----------------------------|
| 1 | Raspas KAREY KAREY Rasps |
| Ref. F0005209 | Número 9 Number 9 |
| Ref. F0005210 | Número 10 Number 10 |
| Ref. F0005211 | Número 11 Number 11 |
| Ref. F0005212 | Número 12 Number 12 |
| Ref. F0005213 | Número 13 Number 13 |
| Ref. F0005214 | Número 14 Number 14 |
| Ref. F0005215 | Número 15 Number 15 |

| | |
|---------------|---|
| 2 | Cabezas femorales de prueba para rasps Femoral trial heads for rasps |
| Ref. A1550040 | Ø 28 C/C Short Neck |
| Ref. A1550041 | Ø 28 C/M Medium Neck |
| Ref. A1550042 | Ø 28 C/L Long Neck |
| Ref. A1550043 | Ø 28 C/XL Extra Long Neck |
| Ref. A1550050 | Ø 26 C/C Short Neck |
| Ref. A1550051 | Ø 26 C/M Medium Neck |
| Ref. A1550052 | Ø 26 C/L Long Neck |
| Ref. A1550053 | Ø 26 C/XL Extra Long Neck |

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| | | | | | |
|---------------|----------------------|---------------|---------------------------|---------------|---------------------------|
| Ref. A1536040 | Ø 28 C/C Short Neck | Ref. A1536043 | Ø 28 C/XL Extra Long Neck | Ref. A1536052 | Ø 26 C/L Long Neck |
| Ref. A1536041 | Ø 28 C/M Medium Neck | Ref. A1536050 | Ø 26 C/C Short Neck | Ref. A1536053 | Ø 26 C/XL Extra Long Neck |
| Ref. A1536042 | Ø 28 C/L Long Neck | Ref. A1536051 | Ø 26 C/M Medium Neck | | |

Bandeja inferior
Lower tray

| | |
|---------------|---|
| 1 | Impactor reductor Impactor-reducer |
| Ref. A2700375 | |
| 2 | Fresa iniciadora endomedular Endomedullary starting reamer |
| Ref. A2700600 | |
| 3 | Cinzel hueco nº 2 Hollow box chisel nº 2 |
| Ref. A2300805 | |
| 4 | Extractor para cabeza femoral Femoral head extractor |
| Ref. A1150010 | |

| | |
|--------------------|--|
| 5 | Punzón impactador para vástago Impactor awl for stem |
| Ref. F0005500 | |
| 6 | Mango impactor vástago definitivo Impactor handle definitive stem |
| Ref. F0005300 | |
| 7 | Varilla para mango de raspa Tommy bar for rasp handle |
| Ref. A1530101 | |
| 8 | Mango automático para rasps Automatic handles for rasps |
| Ref. A1555005 | |
| 2 unidades 2 units | |

Implantes Implants


Karey-HA

| | |
|---------|----------------|
| Num. 9 | Ref. F0005009E |
| Num. 10 | Ref. F0005010E |
| Num. 11 | Ref. F0005011E |
| Num. 12 | Ref. F0005012E |
| Num. 13 | Ref. F0005013E |
| Num. 14 | Ref. F0005014E |
| Num. 15 | Ref. F0005015E |

Karey-C

| | |
|---------|----------------|
| Num. 9 | Ref. F0005109E |
| Num. 10 | Ref. F0005110E |
| Num. 11 | Ref. F0005111E |
| Num. 12 | Ref. F0005112E |
| Num. 13 | Ref. F0005113E |
| Num. 14 | Ref. F0005114E |
| Num. 15 | Ref. F0005115E |

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Cabezas femorales de acero inoxidable*
Stainless-steel femoral heads*

| | |
|----------------|---------------------------|
| Ref. A1507040E | Ø 28 C/C Short Neck |
| Ref. A1507041E | Ø 28 C/M Medium Neck |
| Ref. A1507042E | Ø 28 C/L Long Neck |
| Ref. A1507050E | Ø 26 C/C Short Neck |
| Ref. A1507051E | Ø 26 C/M Medium Neck |
| Ref. A1507052E | Ø 26 C/L Long Neck |
| Ref. A1509040E | Ø 28 C/C Short Neck |
| Ref. A1509041E | Ø 28 C/M Medium Neck |
| Ref. A1509042E | Ø 28 C/L Long Neck |
| Ref. A1509043E | Ø 28 C/XL Extra Long Neck |
| Ref. A1509050E | Ø 26 C/C Short Neck |
| Ref. A1509051E | Ø 26 C/M Medium Neck |
| Ref. A1509052E | Ø 26 C/L Long Neck |
| Ref. A1509053E | Ø 26 C/XL Extra Long Neck |

*Disponible cabeza femoral en Cromo Cobalto
*Availble femoral head Cobalt-Chrome

Distribuidor oficial de

Tecnología Europea de Vanguardia
Advanced European Technology

Industria del Agave 159-A
Zapopan Industrial Norte
45130, Zapopan, Jalisco
33 2471 0012

spinetech@spinetech.com.mx

www.spinetech.com.mx / www.surgival.com